



Dragon Medical Practice Edition, Version 2.3

Release Notes



Copyright

2016. Dragon Medical Practice Edition

Copyright © 2002-2016 Nuance Communications, Inc. All rights reserved.

Nuance, ScanSoft, the Nuance logo, the Dragon logo, Dragon, DragonBar, NaturallySpeaking, NaturallyMobile, RealSpeak, Nothing But Speech (NBS), Natural Language Technology, Select-and-Say, MouseGrid, and Vocabulary Editor are registered trademarks or trademarks of Nuance Communications, Inc. in the United States or other countries. All other names and trademarks referenced herein are trademarks of Nuance Communications or their respective owners. Designations used by third-party manufacturers and sellers to distinguish their products may be claimed as trademarks by those third-parties.

Updated 6/27/2016

L-3892

Contents

- Overview of Dragon Medical Practice Edition, Version 2.3 2
- About this release..... 2
- Upgrading Dragon Medical Practice Edition..... 3
- What's new in Dragon Medical Practice Edition 2.3 4
- Known issues and fixes in Dragon Medical Practice Edition 2.3..... 6
- Product build numbers 7

Overview of Dragon Medical Practice Edition, Version 2.3

Version 2.3 is the newest US English release of Dragon Medical Practice Edition, Nuance's leading speech-recognition and reporting solution for medical professionals. It is released for US English only, and provides enhancements and fixes for current customers of the US English product.

Version 2.3 is available only as a service pack upgrade, downloaded from the Dragon Software Manager.

About this release

Product:	Dragon Medical Practice Edition, version 2.3
Software version:	12.53.350.033 US English, service pack For a comprehensive list of Dragon Medical Practice Edition software versions, see " Product build numbers " on page 7.
Where to download:	Download the update from the Dragon Software Manager. For instructions, see " Upgrading Dragon Medical Practice Edition " on page 3.
Software components provided:	This release is made available via the Dragon Medical Practice Edition 2.3 Service Pack installer, DMPE2DOT3.exe . When installation is complete, these additional components are copied to folders in the program directory: <ul style="list-style-type: none"> • DragonClientForRemoteDesktop • Citrix extensions: <ul style="list-style-type: none"> – vddnspatch2.exe – vddnspatch2.msi – PowerMic Citrix Client Extension.msi – PowerMic Citrix Server & Virtual Desktop Extension.exe Use these files, as appropriate, to install the components you require.
Supported upgrade paths:	Customers currently running version 2.2 (build 12.52.350.048) of Dragon Medical Practice Edition can upgrade to version 2.3 with this service pack. There is no direct upgrade path from Dragon Medical Practice Edition version 1.0 (11.00.304.530 or 11.00.300.312) or from Dragon Medical Practice Edition 2.x versions older than version 2.2. If you are using a version older than version 2.2, you must first upgrade to Dragon Medical Practice Edition 2.2 before you can install the Dragon Medical Practice Edition 2.3 service pack.
Release-specific comments and notes:	Dragon Medical Practice Edition 2.3 adds support for: <ul style="list-style-type: none"> • Windows 10 • Microsoft Office 2016 • Nuance PowerMic III • New edit controls • Internet Explorer 11 • Citrix XenApp 7.6 and 7.7 • Citrix XenDesktop 7.6 and 7.7 • Latest Citrix Receiver 4.4.x
Software documentation:	To view the latest documentation, see the Dragon Medical Practice edition guides in solution 14526 on iSupport .

Upgrading Dragon Medical Practice Edition

The version 2.3 upgrade is available through the Dragon Software Manager.

Notes:

- This service pack only supports upgrading from Dragon Medical Practice Edition version 2.2.
- This service pack does not include updates for language or acoustic models.

About the Dragon Software Manager

Depending on your settings, the Dragon Software Manager may automatically alert you that an update is available for you to install.

If you receive the alert while you are in the middle of a task to modify your Dragon installation—such as creating a new profile, adding words to the vocabulary, or configuring formatting options—please complete this work and save your changes before installing the update. Otherwise, files the update needs to access may be locked and the update would fail.

Installing the update

To upgrade to version 2.3 of Dragon Medical Practice Edition, follow these steps:

1. Launch your existing Dragon Medical Practice Edition software and select **DragonBar>Help>Check for Updates** to open the Dragon Software Manager.
2. (Optional) To display the content of this release notes document on your screen before proceeding, click the line for the Dragon Medical Practice Edition 2.3 Service Pack update to expand the entry, then click the **Read More** link.
3. Select the checkbox for the 2.3 update in the "New Updates & Messages" field of the Dragon Software Manager, and do one of the following:
 - To install the upgrade now, click **Install**.
 - To download the service pack so you can install it later, click **Download Only**. The Dragon Software Manager saves the file as DMPE2DOT3.exe. Double-click it when you're ready to install the upgrade.

Additional software components

The 2.3 upgrade installer (DMPE2DOT3.exe) does not automatically install the Citrix Extensions or the Dragon Client for Remote Desktop updates, although it does copy the files you need for installing these components into the same directory where the main Dragon executable is installed.

If you need to install these components, you can find the files you need within the appropriate subdirectories. For example, in Windows 8 (64-bit), the default location for these components is:

`C:\Program Files (x86)\Nuance\NaturallySpeaking12\Program\Citrix Extension`

and

`C:\Program Files (x86)\Nuance\NaturallySpeaking12\Program\DragonClientForRemoteDesktop`

Note: In the 32-bit version of Windows 8, you'll find these both of these subdirectories in

`C:\Program Files\Nuance\NaturallySpeaking12\Program\`

What's new in Dragon Medical Practice Edition 2.3

Dragon Medical Practice Edition 2.3 adds functionality and support for the latest hardware, software, and operating systems, as summarized in the table below.



New in this release
<p>Windows 10 Dragon Medical Practice Edition 2.3 supports the Windows 10 operating system.</p>
<p>Microsoft Office 2016 Dragon Medical Practice Edition 2.3 supports Microsoft Office 2016.</p>
<p>Nuance PowerMic III Dragon Medical Practice Edition 2.3 supports the PowerMic III. For more information, see "Using the Nuance PowerMic III" on page 4.</p>
<p>Internet Explorer 11 Dragon Medical Practice Edition 2.3 supports Internet Explorer 11.</p>
<p>Support for edit controls (Administrators only) Dragon Medical Practice Edition 2.3 provides support for the following edit controls:</p> <ul style="list-style-type: none"> • TER22 • TER23 • EpicTer83Class


Using the Nuance PowerMic III

Dragon Medical Practice Edition 2.3 supports the Nuance PowerMic III, including new functionality for configuring its buttons, which Dragon enables automatically when you plug in the device.

To view or change the button settings, do the following:

1. Select **DragonBar>Tools>Options** to open the Options dialog, and click the PowerMic tab.
2. Click **Additional** to open the Additional PowerMic Options dialog.
3. Configure the buttons to suit your needs. By default, all buttons are set to behave like a mouse or touchpad. You can choose, instead, to have Dragon execute For each button, you can:
 - Assign it to execute a command. If necessary, clear (unselect) the **Mouse Action** check box for the button, then select Global and Dictation Box commands that you want Dragon to execute when you press the button.
 - Enable Mouse Action. Do this if you want to use the button like a mouse button. Each button has just one mouse action, described in the table below. You can disable a button's mouse action but you cannot change it. When you enable mouse action for a button, you can't also use that button to execute a command, and you can't enable mouse action only for the Dictation Box.

Setting name	Button	Behavior
Mouse Left or Touchpad press		Same as clicking left mouse button. Both buttons are enabled with this check box.
Mouse Right		Same as right-clicking with a mouse; usually will open context menu.

Trigger (underneath) <i>This button is on the back of the PowerMic.</i>		Same as clicking left mouse button.
--	---	-------------------------------------

4. (Optional) Select **Enable Touchpad swipe** if you want to use the PowerMic's touchpad button as a touchpad that you can swipe to move the mouse pointer.
5. Click **OK** to save your changes or click **Restore defaults** to reset the screen to the default button configuration.

Known issues and fixes in Dragon Medical Practice Edition 2.3

The following issue was fixed in Dragon Medical Practice Edition, version 2.3.

Issue	In-Depth Description
Issue 1108 - DragonBar and Mic freeze	<p>Issue Too many internal messages overwhelm the system and cause Dragon to appear frozen.</p> <p>Solution Optimized message handling to avoid lockup. DragonBar and Mic no longer appear frozen.</p>

Product build numbers

The following tables summarize the builds that are part of different releases of Dragon Medical Practice Edition.

Dragon Medical Practice Edition, Version 2.x Releases

Product version	Build Number
2.3 Available as a service pack only.	12.53.350.033
2.2 Available as both a full build and as a service pack.	12.52.350.048
2.1 Available as a service pack only.	12.51.350.006
2.0 Available as a full build only.	12.50.350.022

Dragon Medical Practice Edition, Version 1.0 Releases

Build Number
11.00.304.530
11.00.300.312

How to find the Dragon Medical Practice Edition build number

Method 1:

On the DragonBar, click **Help>About Dragon**.

Method 2:

Use Windows Explorer to locate and open the dragon.log file, where you will find the build number beside the text "Dragon Medical Version".

- In Windows Vista, Windows 7, Windows 8.x, Windows 10, Windows Server 2008, or Windows Server 2012 and 2012 R2 the file is located at:

C:\Users\<<WindowsUsername>\AppData\Roaming\Nuance\NaturallySpeaking12\Dragon.log

- In Windows XP, the file is located at:

C:\Documents and Settings\<<user_account_login>\ApplicationData\Nuance\NaturallySpeaking12\Dragon.log