DRAGON MEDICAL PRACTICE EDITION

QUICK REFERENCE CARD




Dragon Medical Practice Edition 2

Quick Reference Card

CONNECT YOUR MICROPHONE – important notes:

• When you plug your microphone into your PC, an "audio event" window may open. If this happens, verify what is highlighted in that window before closing it.

• If you want to bypass your PC's sound system (or if your PC has only one audio socket), consider using a **USB adapter** or **USB microphone**.



• If your PC has only one audio socket, you can use a Y-shaped adapter.



- If you use a Bluetooth microphone, plug its dongle into a USB port on your PC. After powering on the headset, make sure it is in call mode (see your microphone's guide).
- Using a **built-in (array)** microphone with Dragon Medical Practice Edition is not recommended.
- You may want to verify the "recording" Sound settings in the **Control Panel** (right-click the speaker icon in the System Tray).



IMPORTANT: If your microphone requires plugging into a USB port, aim to always use the same port and plug in before you open Dragon Medical Practice Edition.

INSTALL DRAGON MEDICAL PRACTICE EDITION

NOTE: If you have an existing version, we suggest backing up your profile first.

- 1) Exit any open program. If you're running an anti-virus program, turn it off until after the installation completes.
- 2) Put the Dragon DVD into your PC (or open the downloaded file) and follow the prompts.

3) Enter the **serial number** found on the DVD sleeve (or, if you purchased a download, in the email you received). Click Next, accepting all the recommended options.

Once installation is finished, make sure your microphone is connected, then startDragon Medical Practice Edition (double-click its **desktop icon**).



CREATE OR UPGRADE A USER PROFILE

IF YOU ARE NEW TO DRAGON MEDICAL PRACTICE EDITION... A series of screens will guide you to create the best possible profile for you. Don't rush through it. It will determine your initial accuracy!

Once your profile is created, Dragon Medical Practice Edition offers ways to refine it, including designating **specific text to analyze**.

If you used Dragon Medical 10, 11 or Dragon Medical Practice Edition 1...

Dragon Medical Practice Edition will detect profiles on your PC and start its Upgrade Wizard. (If you do not want to upgrade any profile now, close the wizard. You can launch it later from the Start menu.)

IMPORTANT: Once the upgrade completes, we strongly recommend you **launch Accuracy Tuning** (see the **Audio** menu) and read the **What's New**.



BEFORE YOU START DICTATING...

Go through the **Interactive Tutorial** (its simulations include practice for efficient editing including the **"go back"** command and **smart formatting**).

Then click the DragonBar's Help menu and see its content, especially the regular Tutorial.



PERSONALIZE YOUR PROFILE'S VOCABULARY

In order to dictate efficiently, you should let Dragon Medical Practice Edition learn what you use: nicknames, jargon, capitalized phrases, email addresses...

Personalizing your Profile will prevent many errors.

There are several means to do it including having Dragon Medical Practice Edition analyze text similar to what you plan to dictate. See "Improve my accuracy" in the Help menu.

:

QUICK REFERENCE CARD, CONT.

Dragon Sidebar X

Global

listening"
"Wake up" or "Listen to me"
"Microphone off"

"What can I say?"
"Dragon Sidebar Home"
"Dragon Sidebar Mouse"
"Dragon Sidebar dock right"

"Dragon Sidebar float" "Dragon Sidebar wider

Global Commands

THE DRAGON SIDEBAR

The Sidebar displays tips and sample commands (global and application-specific).

The top pane contains clickable links to related Help topics and tools

You can open the Sidebar by saying "Open Dragon Sidebar" or "What Can I Say?"



NOTE: Right-clicking the Sidebar displays a menu from which you can hide the Mouse tab, print, and more.

AN IMPORTANT HABIT: CONTROL WHETHER DRAGON LISTENS

One can control Dragon's microphone by:

- Pressing the + on the numeric keypad (this hotkey can be changed in Dragon's Options dialog).
- Clicking the microphone icon on the DragonBar or System Tray.
- Using commands such as "wake up" and "microphone off."
 On: Dragon is listening.
 Off: Dragon cannot hear anything.





YOUR FIRST DICTATION

Start a word processor and make sure your cursor is in the document, then turn on the microphone and dictate sentences in a clear and natural voice. Say punctuation and commands like **new paragraph**.

As you dictate, a small Dragon icon indicates that the software is processing. No need to wait for Dragon Medical Practice Edition to "catch up." Dictate at a natural pace. Aim to speak long phrases, since this gives Dragon Medical Practice Edition some context and leads to much better results than slow, halting speech.

NOTE: Dragon Medical Practice Edition automatically handles a lot of formatting, including spacing, dates and percentages. If you want to capitalize a particular word, say "cap" or "all 3

caps" before it. CORRECTING ERRORS

Although personalizing the vocabulary will prevent many errors, some errors are inevitable. By correcting Dragon Medical Practice Edition's errors, you can help it learn from them. Experiment with *when* and *how* to correct. Dragon Medical Practice Edition lets you do it entirely by voice, but also allows **typing over**, pressing a **hotkey**, and more. (See the Options dialog.)

When Dragon Medical Practice Edition misrecognizes something you dictated, you can say "Correct that" or "Correct <xyz>" -- <xyz> being the wrong word(s) on the screen.

This will display the **Correction Menu**, which shows a list of alternatives. If what you wanted is on the list, choose it (say for instance "choose 1").

If not, you can say "Spell that" to open the Spelling window, where you can **spell or type** what Dragon should have recognized or **modify** an alternative that is almost correct.



SHUTTING DOWN DRAGON

It's useful to **exit Dragon Medical Practice Edition** at least once a day: from the **Profile menu**, choose Exit Dragon.


If a message asks whether you want to save changes to your profile, say yes to preserve the refinements made during your session (such as analyzing text or adding a Spoken Form). Note: Dragon may offer to also perform a **quick optimization** of your profile.

LEARNING MORE

At first, focus on getting used to dictating fluently and **refining the vocabulary**. Over time, you can learn to do more and more by voice.

At any time, you can search the Help with a direct command such as "Search Dragon Help for hotkey options." You can also use Dragon Medical Practice Edition's Help menu.

Many resources are on the Web including the latest guides, tips, videos, compatible hardware, partners who provide customization services and Support.



4

QUICK REFERENCE CARD, CONT.

SAMPLE CHARACTERS (see the Help for more Medical Practice Edition) NOTE: You can adjust spoken forms and properties in the Vocabulary Editor. Vocabulary Editor Spoken form (f different) new Ine spoken form al caps of al caps of al caps off al caps off caps off caps off caps off no papace off Norder-Words with spoken forms only Vou Can Say... To Enter...

dot | point | period | full stop exclamation point | mark apostrophe ess open paren close quote at sign (a) plus sign underscore hyphen semi colon (s) optional ess winky face ;-) & ampersand section sign hashmark | number sign

5

NUMBERS, DATES, TIMES, UNITS, PRICES...

To Enter...
4.2

Boston, MA 02460

1 Wayside Dr.

\$4.50

You Can Say...

four point two

four dollars and fifty cents

NOTE: You can change the way Dragon Medical Practice Edition formats numbers, abbreviations, and more. (The Vocabulary Editor lets you set **Word Properties** such as "Fig." for "figure" in front of numbers and "colour" for "color.")

4	numeral four (or say spell four)
MMIV	Roman numeral two thousand four
0.03	zero point zero three
4 x 4	four by four
4 mm	four millimeters
4°	four degrees
179	one hundred seventy nine ${\it or}$ one seventy nine ${\it or}$ one hundred and seventy-nine
5423	five thousand four hundred twenty three
11/32	eleven over thirty two
781-565-5000	seven eight one five six five five thousand
8:30 PM	eight thirty P M
200 kg	two hundred kilograms
5′3″	five feet three inches
May 15, 2013	May fifteenth two thousand thirteen ("comma" is optional)



Boston Massachusetts oh two four six oh

One Wayside Drive

6

QUICK REFERENCE CARD, CONT.

SAMPLE COMMANDS (see Sidebar, Help and Command Browser)

- To click a button or other interface item, say its **name preceded by "click"** (see the Commands tab of the Options dialog).
- Pause before and after commands but not within them.

Control the Microphone Go to sleep | Stop listening

Wake up Microphone off

Search the Web Open top site for...

Search Google for 7 divided by 12 Search Wikipedia for... Search maps for...

Search news for...

Select Text Select all

Select <xvz>

Select next <n> words Select <start> through <end>

Select previous paragraph Select document

Unselect that

Correct Dragon's Errors

Correct <xvz> Correct that

Move the Insertion Point Insert before <xyz>

Go back Go to top | bottom Move down <n> lines Go to end of line Move left <n> characters

Page up | down

Get Help

Give me help

What can I say Dragon Sidebar Help Search Dragon Help for...

Search the Computer

(requires a desktop search engine) Search the computer for... Search documents for... Search email for...

Edit Text

Resume with <xyz>

Delete line

Delete last <n> words Scratch that <n> times

Backspace <n>

Undo that Cut that Paste that Show Dictation Box

Spell Out Spell that Spell <cap b a hyphen 5>

Spell <space Charlie alpha papa> Switch to Spell mode

Add Lines and Spaces

New line New paragraph Press Enter Press Tab key Tab <n> times

Move in a list Format

Move down <n> Bullet selection, Unbullet that Go to bottom | top Bold the previous line

Underline <xyz>, Capitalize <xyz> Press Enter Press right arrow Make this uppercase, Uncap that

All caps on | off

Work with Windows

Switch to <window name> Minimize window

Show Desktop Restore windows List all windows List windows for cprogram>

Start and Close Items (see the options for Start menu and Desktop)

Start <item name> Start Microsoft Word Start Mail Start Internet Explorer Close window, Press alt F4 Open Control Panel

Move the Mouse Position the Mouse Move mouse up MouseGrid Move mouse right MouseGrid window MouseGrid <1 to 9><1 to 9> Move mouse down slower

Click the Mouse Drag the Mouse Mouse click Drag mouse down faster

Mouse double click Mouse drag lower right very fast Mouse right click Mouse drag up very fast

Medical templates Medical normals Dragon brief produce template Dragon brief female exam Dragon inpatient admit template Dragon brief male exam Dragon new patient template Dragon normal chest Dragon progress note template Dragon normal heart

Dragon discharge summary template Dragon normal skin



Copyright © 2013 Nuance Communications, Inc. All rights reserved. Nuance, Dragon, the Dragon $logo, and \ the \ Dragon Bar \ are \ trademarks \ or \ registered \ trademarks \ of \ Nuance \ Communications, \ Inc. \ or \ Nuance \ Communications \ Property \ P$ its affiliates in the United States and/or other countries. All other names and trademarks referenced NUANCE herein are the trademarks or registered trademarks of their respective owners.

51-A71A-28000