

# ODMS R7.2 Patch Program Release Note

OLYMPUS CORPORATION  
Image Processing Technology 2

May 28<sup>th</sup>, 2019

## ■ 1 ODMS R7.2 Patch Program

The newest patch includes the updates from the previous patch releases listed below.

Patch No.	Description
Patch 0004	<ol style="list-style-type: none"><li>1. Fixed the issue that ODMS cannot find Dragon profile in some cases.</li><li>2. Fixed the issue that ODMS TM installed in Workgroup Mode cannot open the Dragon Profile when executing voice recognition.</li></ol> <p>Remark: After applying this patch program, you need to setup the Dragon profile in ODMS again and restart ODMS in order to have the SCP resynchronize the profile.</p>
Patch 0003	<ol style="list-style-type: none"><li>1. Fixed the bug in Patch0002</li></ol>
Patch 0002	<ol style="list-style-type: none"><li>1. Fixed the issue that USB devices such as RecMic, DVR or Foot Switch may not be recognized by the application during startup.</li><li>2. Fixed the issue that USB devices such as RecMic, DVR or Foot Switch cannot be recognized in some environment.</li></ol>
Patch 0001	<ol style="list-style-type: none"><li>1. Fixed the issue that if you change the order of the header column by drag &amp; drop in the list view, the change will not be saved in ODMS.</li><li>2. Fixed the icon status error for the In-line Edit function.</li><li>3. Fixed the issue for the device settings for DS-7000 and DS-3500.</li><li>4. Fixed the issue that, when you use the older generation of the Foot Switch (RS27, RS28) by USB connection, it takes time to complete the initialization. There had been cases in a Virtual Environment where you cannot login since the confirmation of connection had been taking time.</li></ol>

[EOF]